

Softclix Lancet + Accu-chek Guide meter



Set-up:

Turn the cap to the number you want.
1 is for sensitive skin
5 is for thicker skin
2-3-4 are inbetween



1. Pull cap straight off



2. Insert lancet...
then twist off round
cover



3. Line up notches
and replace cap.
You will hear a click



4. Press the button all
the way in (like a
ballpoint pen)



5. Touch side of
finger to end of
lancing device




6. Press yellow
button to prick your
finger. You will hear a
click



7. Gently squeeze
finger to get a drop of
blood



8. Remove one test strip



9. Insert gold side into meter.



10. Push strip all the way in. Meter will turn on automatically



11. Apply blood at the yellow edge of the test strip



12. See your result. Pull test strip out to turn meter off



13. Remove cap



14. Use body of lancing device to push lancet into sharps box. Replace cap



Done!