



InPen has a one-year warranty and does not require charging.

## QUICK START GUIDE

Smart Insulin Pen + App system that gives you the tools to manage your multiple daily injections.



**Your personal therapy settings are needed before first use of the InPen dose calculator.**

For additional information, refer to the Instructions for Use in the InPen App.

**STEPS TO BEGIN USING YOUR INPEN**

- 1** **Download the InPen App and Create an Account:**  
Open App Store and search for InPen. Download InPen by Companion Medical App. Create login with email, Facebook, Gmail, or Twitter.



- 2** **Sign in and accept the Terms of Service.**

- 3** **Follow the setup wizard to pair the InPen to your smartphone.**

- 4** **Need support? We are here to help.** Contact Customer Support at 844.843.7903 or support@companionmedical.com.

**1 YEAR WARRANTY**

Join the community



To learn more visit: [www.companionmedical.com](http://www.companionmedical.com)  
or call 844.843.7903